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6. Yamaguchi T, Nakamura I, Chiba K, et al. Epidemiological and microbiological analysis of community-associated methicillin-resistant *Staphylococcus aureus* strains isolated from a Japanese hospital. *Jpn J Infect Dis.* 65:175-178. 2012.
7. Inomata S, Yano H, Tokuda K, et al. Microbiological and molecular epidemiological analyses of community-associated methicillin-resistant *Staphylococcus aureus* at a tertiary care hospital in Japan. *J Infect Chemother.* 21:729-736. 2015.
8. Paterson GK, Harrison EM, Holmes MA. The emergence of *mecC* methicillin-resistant *Staphylococcus aureus*. *Trends Microbiol.* 42:7. 2014.
9. Bal AM, Coombs GW, Holden MTG, et al. Genomic insights into the emergence and spread of international clones of healthcare-, community- and livestock-associated methicillin-resistant *Staphylococcus aureus*: Blurring of the traditional definitions. *J Glob Antimicrob Resist.* 6:95-101. 2016.
10. Clinical and Laboratory Standards Institute (CLSI) 2023. Performance standards for antimicrobial susceptibility testing; M100-Ed33.
11. 松本竹久. 臨床検査で問題になる Small-colony variants. *日本臨床微生物学会雑誌*, 33(2), 101-108, 2023.
12. Kipp F, Becker K, Peters G, et al. Evaluation of different methods to detect methicillin resistance in small-colony variants of *Staphylococcus aureus*. *J Clin Microbiol.* 42(3):1277-9. 2004.
13. Clinical and Laboratory Standards Institute (CLSI). Performance Standards for Antimicrobial Disk Susceptibility Tests; M02, 13th Edition. 2015.
14. 国立感染症研究所. 病原体検出マニュアル. 侵襲性インフルエンザ菌感染症 (2021年12月版) . https://www.niid.go.jp/niid/images/lab-manual/Invasive_Haemophilus_influenzae_disease20211228.pdf (参照 2024-02-22)
15. Clinical and Laboratory Standards Institute (CLSI). Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; M07, 11th Edition. 2018.

16. 石井良和. 薬剤感受性試験とブレイクポイント, その問題点と今後の展望. 日本化学療法学会雑誌, 59, 454-459, 2011.
17. Kato F, Nakamura M, Sugai M. The development of fluorescent protein tracing vectors for multicolor imaging of clinically isolated *Staphylococcus aureus*. *Sci Rep.* 7:2865. 2017.
18. Petersdorf S, Herma M, Rosenblatt M, et al. A Novel Staphylococcal Cassette Chromosome mec Type XI Primer for Detection of mecC-Harboring Methicillin-Resistant *Staphylococcus aureus* Directly from Screening Specimens. *J Clin Microbiol.* 53(12):3938-3941. 2015.
19. 薬剤耐性菌バンク. <https://jarbb.jp> (参照 2024-02-22)
20. Kondo Y, Ito T, Ma XX, et al. Combination of multiplex PCRs for staphylococcal cassette chromosome mec type assignment: rapid identification system for mec, ccr, and major differences in junkyard regions. *Antimicrob Agents Chemother.* 51:264-274. 2007.
21. Zhang K, Sparling J, Chow BL, et al. New quadriplex PCR assay for detection of methicillin and mupirocin resistance and simultaneous discrimination of *Staphylococcus aureus* from coagulase-negative staphylococci. *J Clin Microbiol* 42:4947-4955. 2004.
22. Kawamura H, Tokuda K, Imuta N, et al. Epidemiological Analysis of Nosocomial MRSA Outbreaks Using Phage Open-Reading Frame Typing in a Tertiary-Care Hospital. *Jpn J Infect Dis.* 69:523-524. 2016.

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