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8. Inomata S, Yano H, Tokuda K, et al. Microbiological and molecular epidemiological analyses of community-associated methicillin-resistant *Staphylococcus aureus* at a tertiary care hospital in Japan. *J Infect Chemother.* 21:729-736. 2015.
9. CLSI. *Performance Standards for Antimicrobial Susceptibility Testing.* 36th ed. CLSI supplement M100. Clinical and Laboratory Standards Institute; 2026.
10. 松本竹久. 臨床検査で問題になる Small-colony variants. *日本臨床微生物学会雑誌*, 33(2), 101-108, 2023.
11. Kipp F, Becker K, Peters G, et al. Evaluation of different methods to detect methicillin resistance in small-colony variants of *Staphylococcus aureus*. *J Clin Microbiol.* 42(3):1277-9. 2004.
12. CLSI. *Performance Standards for Antimicrobial Disk Susceptibility Tests.* 14th ed. CLSI standard M02. Clinical and Laboratory Standards Institute; 2024.
13. CLSI. *Performance Standards for Antimicrobial Susceptibility Testing; Twenty-Second Informational Supplement.* CLSI document M100-S22. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
14. CLSI Archived Resource. A reference of breakpoints eliminated from CLSI document CLSI M100 Since 2010.
https://clsi.org/media/bjhnzhh/m100_archived_drugs_table_20251022.pdf
15. 国立感染症研究所. 病原体検出マニュアル. 侵襲性インフルエンザ菌感染症 (2021年12月版). https://id-info.jihs.go.jp/manuals/pathogen-detection/Invasive_Haemophilus_influenzae_disease20211228.pdf (参照 2026-6-9)
16. CLSI. *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically.* 12th ed. CLSI standard M07. Clinical and Laboratory Standards Institute; 2024.
17. 石井良和. 薬剤感受性試験とブレイクポイント, その問題点と今後の展望. *日本化学療法学会雑誌*, 59, 454-459, 2011.

18. Kato F, Nakamura M, Sugai M. The development of fluorescent protein tracing vectors for multicolor imaging of clinically isolated *Staphylococcus aureus*. Sci Rep. 7:2865. 2017.
19. 薬剤耐性菌バンク. <https://jarbb.jp> (参照 2024-02-22)
20. Kondo Y, Ito T, Ma XX, et al. Combination of multiplex PCRs for staphylococcal cassette chromosome *mec* type assignment: rapid identification system for *mec*, *ccr*, and major differences in junkyard regions. Antimicrob Agents Chemother. 51:264-274. 2007.
21. Petersdorf S, Herma M, Rosenblatt M, et al. A Novel Staphylococcal Cassette Chromosome *mec* Type XI Primer for Detection of *mecC*-Harboring Methicillin-Resistant *Staphylococcus aureus* Directly from Screening Specimens. J Clin Microbiol. 53(12):3938-3941. 2015.
22. Zhang K, Sparling J, Chow BL, et al. New quadriplex PCR assay for detection of methicillin and mupirocin resistance and simultaneous discrimination of *Staphylococcus aureus* from coagulase-negative staphylococci. J Clin Microbiol 42:4947–4955. 2004.
23. Kawamura H, Tokuda K, Imuta N, et al. Epidemiological Analysis of Nosocomial MRSA Outbreaks Using Phage Open-Reading Frame Typing in a Tertiary-Care Hospital. Jpn J Infect Dis. 69:523-524. 2016.
24. Kaya H, et al. SCCmecFinder, a Web-Based Tool for Typing of Staphylococcal Cassette Chromosome *mec* in *Staphylococcus aureus* Using Whole-Genome Sequence Data. mSphere. 3:e00612-17. 2018.
25. Petit RA III. sccmec: A tool for typing SCCmec cassettes in assemblies. <https://github.com/rpetit3/sccmec>

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